NEWTON 7.0 FT-500

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CAMERA & OPTICS

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

Applications 应用实例

Bioluminescence-肿瘤3D成像

4T1 mouse breast cancer cell line carrying *Luciola Italia* (Red-Fluc) luciferase, Image was captured on the day at 1 week after injection of cancer cells and signal was located in the mammary gland. 2D bioluminescence image and corresponding topography, tomography images show.


Bioluminescence-肿瘤发育研究

Locoregional therapy with astatine-211-labeled trastuzumab ($^{211}\text{At-tras}$) in a mouse model of peritoneal metastasis of gastric cancer (PMGC), and tumor monitoring. (a) Treatment schedule. (b) Representative bioluminescence images of tumor growth in the mouse model treated with PBS (control), unlabeled trastuzumab (Tras), ^{211}At (1 MBq), or $^{211}\text{At-tras}$ (0.1 or 1 MBq). Images were captured on the day before treatment (Day 0) or at 3, 5, or 7 weeks (w) after treatment. The color scale indicates the luminescence intensity per pixel. *Cancer Sci.* 2017; 108: 1648-1656


Bioluminescence-肿瘤发育研究

In vivo antitumor activity of NK-92- $\alpha\text{FR-28BB}\zeta$ cells. A, Schematic representation of the construction of mouse xenograft model of ovarian cancer and treatment procedures. B, Ovarian cancer development was monitored by in vivo bioluminescence imaging. Images taken at day 0 and day 10 of the experiment are shown. C, Quantification summary of bioluminescence signal intensity in each group shown in (B). The data are expressed as the means \pm SEM of 5 mice for all groups. * $P < 0.05$; ** $P < 0.01$; ns, $P \geq 0.05$. D, Kaplan-Meier survival curves of tumor-bearing mice treated with PBS, NK-92 cells, NK-92-EV cells, or NK-92- $\alpha\text{FR-28BB}\zeta$ cells. ** $P < 0.01$, compared with the NK-92- $\alpha\text{FR-28BB}\zeta$ group. *J Immunother.* 2019;42:284-296


Fluorescence-药物代谢

(A) In vivo noninvasive fluorescence imaging DOX distribution in rats. (B) Ex vivo fluorescence imaging of organs at 10 h post-oral-administration. (C) DOX plasma concentration vs. time profiles of oral administration. (D) Quantification of tissue distribution of DOX in rats ($n = 5$). *Eur. J. Pharm. Biopharm.* 2014; 87:197-207

Fluorescence-纳米材料

(a) Pharmacokinetics of BPDC NCs or free DOX over 24 h after intravenous injection to KM mice; (b) biodistribution of BPDC NCs or free DOX at 24 h postinjection to BALB/c mice bearing 4T1 tumors (** $p < 0.01$ between two groups); (c) in vitro CT images of BPDC NCs; (d) CT intensity values of BPDC NC aqueous dispersion at different concentrations; fluorescence images of BALB/c mice bearing 4T1 tumors before and after the administration of (e) BPDC NCs and (f) Ce6. Major organs and tumors were excised from both groups at 24 h post-injection. *Biomater.Sci.*, 2019, DOI:10.1039/C9BM00351G.

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